

SONICARE AIRFLOSS GUIDE




Philips Sonicare AirFloss

An easier way to clean between teeth.

It takes just 60 seconds to clean between all of your teeth with Sonicare AirFloss, and it's easy to use, too. Here's how to make sure you're correctly using your Sonicare AirFloss so it can thoroughly clean between your teeth.

1. Fill the Reservoir

After the handle is charged, attach the nozzle and fill the reservoir with two teaspoons of water or your favorite mouthwash. The reservoir holds enough for two full uses, so you only need to fill it every other use.



2. Prime the Nozzle

For the first use, fill the reservoir; press the power button to turn Sonicare AirFloss on, then press the activation button approximately six times or until spray comes out of the nozzle.

3. AirFloss Placement

Place the tip between teeth and apply gentle pressure. You can begin wherever you like in your mouth.



4. Begin Cleaning

As the tip rests between your teeth, press the activation button. The innovative technology delivers a powerful burst of air and water that effectively yet gently cleans deep between your teeth.

5. Repeat for Each Interdental Area

Glide the guidance tip along the gum line until you feel it settle between the next set of teeth. Press the activation button again, then repeat as you glide the nozzle.



6. Back Teeth

Don't forget the hard-to-reach spaces between your back teeth. The slim, angled nozzle makes it easier to reach.

7. Turn the Power Off

When you are done, press the power button to turn off Sonicare AirFloss. Don't worry if you forget—the Sonicare AirFloss will automatically turn off after four minutes.

8. Cleaning

After use, remove the nozzle and rinse it with water. If you like, you can use a cotton swab to wipe down the reservoir. We don't recommend using anything other than water or mouthwash inside the device.

Visit Sonicare.com to learn about the benefits of Sonicare.